

ZYNQ ADMIN PANEL

Implementation & Internal Validation Specification

Version 1.3 - Concrete Delivery Requirements

Purpose	Define exactly what the development team must implement, test internally and evidence before the next Admin Panel delivery.
Development team responsibility	Implementation, internal QA, documented PASS/FAIL results and evidence.
ZYNQ / Anna responsibility	Independent acceptance validation after the team has completed and documented its own QA.
Current focus	Operational readiness of the Admin Panel and related Clinic/Doctor panels.
Date	14 July 2026

This document replaces the broad audit format with concrete testable expectations.

1. Purpose and Working Method

Following the team's feedback, the previous readiness audit has been converted into a concrete implementation and internal validation specification. Each item below states the exact workflow to test, the expected result and the evidence required from the development team.

Required delivery sequence

- **1. Development team:** Implement any missing functionality.
- **2. Development team QA:** Run every applicable test in this document before delivery.
- **3. Development team report:** Mark each test PASS, FAIL or NOT APPLICABLE and attach evidence.
- **4. ZYNQ / Anna:** Perform an independent acceptance check on representative workflows after internal QA is complete.
- **5. Closure:** A function is considered complete only when the expected result is demonstrated in the stated environment.

Response rules

- Please do not answer only with "fixed" or "working".
- For each PASS, include the tested environment/build and evidence.
- For each FAIL, include cause, responsible person and target date.
- If a test does not apply to the agreed pilot scope, mark NOT APPLICABLE and explain why.

2. Closed, Excluded and ZYNQ-Owned Items

Item	Current status	Treatment in this document
Global Search core logic	Previously validated by Anna	Not reopened; regression only where admin changes affect it.
Primary Tag inclusion/ranking	Validated	Not reopened.
Admin tag/concept search	Implemented and shown by video; Anna verified Admin-to-App mapping	Not an open delivery item.
AI tag generation for new Treatments	Validated	Not reopened; regression check only.
Translation Management	Delivered and largely validated	Only admin save/publish regression tests remain.
Products	Outside current ZYNQ scope	NOT APPLICABLE. Must not be treated as a blocker.
Completeness of Device-to-Treatment data	ZYNQ data-maturation responsibility	Not a development failure. Team must only prove that saved relationships function correctly.

3. Test Environment and Evidence

ADM-ENV-01 - Identify the test environment

Applies to	React Admin Panel, Clinic Panel, Doctor Panel
------------	---

Team test steps	State the exact URL/build/branch used for internal QA. Confirm whether it is development, staging or production-like.
Expected result	Every reported result can be traced to a named build and environment. No test result is accepted without an environment reference.
Required evidence	Screenshot of environment/build identifier and written build/branch/commit reference.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-ENV-02 - Real data versus mock data

Applies to	All in-scope modules
Team test steps	For each module, confirm whether the displayed data and actions use real backend APIs, seeded staging data, mock data or partially implemented endpoints.
Expected result	The status of every in-scope module is explicit. No mock or local-only flow is represented as production-ready.
Required evidence	Short module matrix: real / seeded / mock / partial / excluded.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-ENV-03 - Role and access control

Applies to	Admin, clinic user, doctor/expert user
Team test steps	Log in with each relevant role. Attempt to open and change functions outside that role.
Expected result	Each role sees only authorised modules and can perform only authorised actions. Direct URLs do not bypass permissions.
Required evidence	Screen recording for one allowed and one blocked action per role.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-ENV-04 - Session and authentication behaviour

Applies to	All panels
Team test steps	Test logout, expired session, refresh after session expiry and direct URL after logout.
Expected result	The user is returned to login safely. Unsaved or protected data is not exposed. Clear messaging is shown.
Required evidence	Screen recording of session-expiry and unauthorised-access behaviour.

Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.
----------------------	--

4. Navigation and List-State Persistence

This is the primary operational workflow. The same pattern must work consistently across all launch-critical list/detail modules.

ADM-NAV-01 - Reference list-to-detail-to-list workflow

Applies to	Customers / Users
Team test steps	Open the customer list. Go to page 3. Apply a search and any available filter/sort. Scroll to a record (reference: record no. 14), open it, make no change, then use Back to list.
Expected result	The exact previous page, page size, search, filters, sort, scroll position and selected-row context are restored.
Required evidence	One continuous screen recording from list setup through return.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-NAV-02 - Cross-module state persistence

Applies to	Clinics, Doctors/Experts, Appointments, Treatments/Services, Devices, Concerns
Team test steps	Repeat ADM-NAV-01 in every applicable module.
Expected result	Each module preserves the same list context. Behaviour is consistent, not module-specific or accidental.
Required evidence	One recording per module or one combined recording clearly showing each module.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-NAV-03 - Browser Back and explicit Back to list

Applies to	All list/detail modules
Team test steps	Open a detail record and test both the browser Back button and the panel's own Back to list control.
Expected result	Both routes return to the correct list context without resetting state or navigating to an unrelated page.
Required evidence	Screen recording of both navigation methods.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-NAV-04 - Scroll and selected-row restoration

Applies to	Long lists
Team test steps	Open a record far down a long list, return to the list and repeat with another record.
Expected result	The list returns to the previous vertical position and the previously opened row is visually identifiable where appropriate.
Required evidence	Screen recording showing the list position before and after.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-NAV-05 - State after save and cancel

Applies to	All editable detail pages
Team test steps	From a filtered/sorted list, open a record. Test Save and return; then repeat with Cancel/Back without saving.
Expected result	Both workflows return to the previous list context. Save persists data; Cancel does not persist changes.
Required evidence	Recording showing both Save and Cancel paths.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

5. Detail Pages, Refresh and Direct Links

ADM-DET-01 - Refresh an existing detail page

Applies to	Customers, Clinics, Doctors, Treatments, Devices, Appointments
Team test steps	Open a detail page and refresh the browser.
Expected result	The same record reloads correctly from its stable ID. No blank page, missing state or crash occurs.
Required evidence	Recording with visible URL before and after refresh.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-DET-02 - Open detail page in a new tab

Applies to	Same modules as ADM-DET-01
Team test steps	Copy the detail URL and open it in a new tab or browser window.
Expected result	The correct record loads without relying on previous route state.

Required evidence	Recording showing copied URL and successful load.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-DET-03 - Invalid or deleted record ID

Applies to	All ID-based detail pages
Team test steps	Open a URL with an invalid ID and a previously deleted test ID.
Expected result	A clear Not found / no longer available state is shown. The page does not crash or expose unrelated data.
Required evidence	Screenshots of both cases.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-DET-04 - Unsaved changes warning

Applies to	Editable detail pages
Team test steps	Edit one or more fields, then navigate away, use Back or close the tab without saving.
Expected result	The user receives a clear warning before losing unsaved changes, or the system uses an agreed safe autosave model.
Required evidence	Recording of the warning and both Stay/Leave choices.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

6. Search, Filters, Sorting and Pagination

ADM-LIST-01 - Search returns the correct records

Applies to	All core list modules
Team test steps	Use one exact name, one partial name and one term with no results.
Expected result	Exact and partial searches return relevant records; no-result searches show a helpful empty state. Search does not require a page reload.
Required evidence	Recording for three searches in at least three launch-critical modules.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-LIST-02 - Filter behaviour

Applies to	Modules with filters
Team test steps	Apply each available filter individually, then combine filters. Clear them afterward.
Expected result	Results match the selected criteria. Clearing filters restores the complete list and a valid page.
Required evidence	Recording showing single, combined and cleared filters.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-LIST-03 - Sorting

Applies to	Core table modules
Team test steps	Sort each supported column ascending and descending. Open a record and return.
Expected result	Rows sort correctly, the direction is visible and the sort remains applied after detail navigation.
Required evidence	Recording showing both directions and return-to-list persistence.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-LIST-04 - Pagination after search or filter

Applies to	Paginated lists
Team test steps	Move to page 3 or later, then apply a restrictive search/filter.
Expected result	Pagination recalculates safely. The user is moved to a valid page; no empty or out-of-range page appears.
Required evidence	Recording showing page before and after filter/search.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-LIST-05 - Page size and total count

Applies to	Paginated lists
Team test steps	Change page size where available and move across several pages.
Expected result	Displayed row count, total count and page controls remain correct. Page size persists during the session and after detail navigation where agreed.
Required evidence	Recording showing page size and total count.

Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.
----------------------	--

ADM-LIST-06 - Rapid input and repeated actions

Applies to	Search fields and table controls
Team test steps	Type quickly, clear and re-enter search text; click filters/sort controls repeatedly at normal user speed.
Expected result	No stale results, flicker loops, duplicate requests causing wrong data or frozen controls. Final UI state matches the final input.
Required evidence	Recording and, if available, browser network/console evidence showing no uncaught errors.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

7. Create, Edit, Save and Form Validation

ADM-CRUD-01 - Create a valid record

Applies to	Customers or test-safe module, Treatments, Devices, Concerns as applicable
Team test steps	Create one new test record using all required fields. Submit once.
Expected result	Exactly one record is created. Success feedback is shown. The record is searchable and opens correctly after refresh.
Required evidence	Recording from form entry through search and refresh.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-CRUD-02 - Required-field validation

Applies to	All create/edit forms
Team test steps	Submit with each required field empty and with invalid formats where relevant.
Expected result	Submission is blocked. Each error appears beside the correct field in clear language. Existing valid data is not lost.
Required evidence	Screenshots or recording of representative validation errors.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-CRUD-03 - Duplicate prevention

Applies to	Treatments, Devices, Concerns, Clinics, Doctors and other uniqueness-sensitive records
Team test steps	Attempt to create a duplicate using the same name/identifier and a case/spacing variation.
Expected result	The system blocks or clearly warns about duplicates according to the agreed rule. It does not silently create conflicting masterdata.
Required evidence	Recording showing duplicate attempt and response.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-CRUD-04 - Edit and persistence

Applies to	All editable modules
Team test steps	Change a text field, status and relationship where applicable. Save, refresh and reopen from the list.
Expected result	All saved changes remain correct after refresh and reopening. Unchanged fields and relationships are not lost.
Required evidence	Recording of edit, save, refresh and reopen.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-CRUD-05 - Cancel does not save

Applies to	All editable modules
Team test steps	Change fields, select Cancel or Back without confirming save, then reopen the record.
Expected result	The original saved values remain. No partial or hidden save occurs.
Required evidence	Recording showing changed values before Cancel and original values after reopen.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-CRUD-06 - Double-submit protection

Applies to	Create/edit forms
Team test steps	Double-click Save/Submit and test repeated clicks during a slow request.
Expected result	Only one request/result is created. The button shows progress and is disabled while processing.
Required evidence	Recording plus network evidence if available.

Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.
----------------------	--

ADM-CRUD-07 - Long text and boundary values

Applies to	Treatments/Services, Clinics, Doctors, Support content and other text-heavy forms
Team test steps	Test maximum reasonable text length, special characters, Swedish characters and line breaks.
Expected result	Content saves and displays without clipping, corruption, broken layout or silent truncation. Any length limit is clearly communicated.
Required evidence	Screenshots before save and after reopen.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

8. Delete, Deactivate, Approval and Operational Safety

ADM-SAFE-01 - Delete confirmation

Applies to	All deletable modules
Team test steps	Choose Delete on a safe test record. Cancel once, then confirm on the second attempt.
Expected result	The confirmation names the record and action. Cancel changes nothing. Confirm deletes only the intended record and shows clear feedback.
Required evidence	Recording of Cancel and Confirm paths.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-SAFE-02 - Dependency protection

Applies to	Treatments, Devices, Concerns, Clinics, Doctors
Team test steps	Attempt to delete a test record that has active relationships.
Expected result	The system blocks unsafe deletion or clearly explains the downstream impact and approved alternative, such as deactivate. No orphan references remain.
Required evidence	Recording showing the dependency message and resulting data integrity.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-SAFE-03 - Deactivate/reactivate

Applies to	Modules supporting status changes
Team test steps	Deactivate an active test record, verify its admin status and user-facing behaviour, then reactivate it.
Expected result	Status changes are explicit, reversible and reflected correctly wherever the record is used.
Required evidence	Recording across Admin and applicable user-facing surface.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-SAFE-04 - Approve/verify/reject workflow

Applies to	Clinics, Doctors/Experts and any approval-based module
Team test steps	Run an approval/verification/rejection workflow on a test record.
Expected result	The UI explains what each action means, asks for confirmation where appropriate, saves the status and records who/when if audit data exists.
Required evidence	Recording of the complete workflow and resulting status.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-SAFE-05 - Publish/unpublish consequence

Applies to	Treatments/Services and other publishable content
Team test steps	Publish a test record, verify visibility, then unpublish or return to draft.
Expected result	The status and consequence are clear before action. User-facing visibility follows the selected status without exposing drafts.
Required evidence	Recording from Admin status change to user-facing verification.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-SAFE-06 - Error during a destructive action

Applies to	Any safe test module
Team test steps	Simulate or test a failed delete/status request in staging.
Expected result	The UI reports failure clearly and leaves the record in its previous valid state. It does not show success when the backend failed.
Required evidence	Recording or controlled test evidence from staging.

Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.
----------------------	--

9. Relationships and Data Integrity

Note: ZYNQ is responsible for completing the actual Device-to-Treatment data mappings. The development team is responsible for proving that the relationship function itself saves, updates and displays correctly.

ADM-REL-01 - Treatment-to-Device relationship save

Applies to	Treatments/Services and Devices
Team test steps	Create or edit a test Treatment and select one or more Devices. Save and reopen both sides.
Expected result	The selected Devices remain connected to the Treatment. Where the reverse view exists, the Treatment is shown under the correct Device.
Required evidence	Recording from selection through reopen and reverse lookup.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-REL-02 - Remove and replace a relationship

Applies to	Treatments/Services and Devices
Team test steps	Remove one Device relationship and add a different Device. Save and refresh.
Expected result	The old relationship is removed, the new one is stored and no unrelated relationship changes.
Required evidence	Recording before and after save/refresh.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-REL-03 - Clinic-to-Doctor relationship

Applies to	Clinics and Doctors/Experts
Team test steps	Add a Doctor to a Clinic, verify both views, then remove the relationship.
Expected result	The relationship is visible and consistent from both relevant views and updates correctly after removal.
Required evidence	Recording of add, verify and remove.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-REL-04 - Clinic-to-Treatment/Service relationship

Applies to	Clinics and Treatments/Services
Team test steps	Assign and unassign a Treatment/Service to a test Clinic.
Expected result	The Clinic offering is stored correctly, appears in the relevant panel/app surface and disappears when unassigned according to the agreed status rules.
Required evidence	Recording across Admin and applicable user-facing surface.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-REL-05 - Relationship preservation during unrelated edits

Applies to	All relational masterdata
Team test steps	Edit an unrelated text field on a record with existing relationships. Save and reopen.
Expected result	Existing relationships remain unchanged unless deliberately edited.
Required evidence	Before/after screenshots or recording.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-REL-06 - Automatic tag/mapping regression

Applies to	New Treatment/Service creation
Team test steps	Create a test Treatment, choose related Device(s), submit and inspect generated tags/mappings using the existing Admin validation function.
Expected result	AI tags are generated, selected Device relationships are saved and the new Treatment is discoverable without manual updates to old entities.
Required evidence	Recording from creation through Admin validation and app/search verification.
Team response	Regression check: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

10. Translation and Publishing Regression

Translation Management has already been delivered. The following tests are regression and operational-safety checks, not a reopening of the completed Global Search work.

ADM-LANG-01 - Save Swedish and English values

Applies to	Treatments/Services, Devices and other translated content
-------------------	---

Team test steps	Enter or edit both language values, save, refresh and reopen.
Expected result	Each value remains in the correct language field. No language overwrites the other.
Required evidence	Recording of edit, save and reopen.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-LANG-02 - AI translation review and edit

Applies to	Translation Management
Team test steps	Generate a translation, edit it before approval and save the reviewed version.
Expected result	The approved edited version is preserved and used; the AI draft does not overwrite it unexpectedly.
Required evidence	Recording of generated text, manual edit and final saved value.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-LANG-03 - Protected terms

Applies to	Devices/brands and protected terminology
Team test steps	Use content containing names such as Candela, Nordlys, Clarity II, Morpheus8 or another configured protected term.
Expected result	Protected names remain unchanged in translation and publication.
Required evidence	Before/after screenshots of both languages.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-LANG-04 - Incomplete translation handling

Applies to	Launch-critical translated content
Team test steps	Leave one required language empty and attempt to publish.
Expected result	The system follows the agreed rule: block publication or show a clear, actionable warning. The user cannot unknowingly publish incomplete content.
Required evidence	Recording of the incomplete-language action and system response.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-LANG-05 - Admin-to-app language display

Applies to	Translated Treatments/Services
Team test steps	Publish an approved test item and view it in Swedish and English app modes.
Expected result	The correct title, description and labels display in each language without unintended mixing when both translations exist.
Required evidence	Side-by-side screenshots or recording of both app languages.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-LANG-06 - Read More / Show More regression

Applies to	Search and long translated content
Team test steps	Open long result content and expanded result lists in both languages.
Expected result	Read More/Read Less and Show More work consistently without clipping, overlap or unexpected reset.
Required evidence	Recording in Swedish and English.
Team response	Regression check: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

11. Loading, Empty, Error and Recovery States**ADM-ERR-01 - Loading state**

Applies to	All API-backed lists and detail pages
Team test steps	Open lists/details under normal and throttled staging conditions.
Expected result	A clear loading state is shown. Controls do not appear active before data is ready. No blank frozen screen is displayed.
Required evidence	Recording under normal and throttled conditions.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-ERR-02 - Empty state

Applies to	All lists
Team test steps	Use filters/search that return zero rows and open a genuinely empty test module if available.
Expected result	The empty state explains that no records match and provides a relevant next action such as clear filters or create new.

Required evidence	Screenshots of filtered-empty and genuinely-empty states.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-ERR-03 - API failure

Applies to	Representative list, detail and save action
Team test steps	Test a controlled backend/network failure in staging.
Expected result	A clear error is shown with Retry where appropriate. Existing data is not silently lost. No false success message appears.
Required evidence	Recording of failure and successful retry.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-ERR-04 - Slow save and navigation lock

Applies to	Representative create/edit form
Team test steps	Submit under slow network conditions and attempt to navigate or submit again.
Expected result	The UI shows progress, prevents duplicate submission and warns before abandoning an in-progress operation if needed.
Required evidence	Recording under throttled conditions.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-ERR-05 - Concurrent edit conflict

Applies to	One record opened in two sessions/tabs
Team test steps	Edit the same record differently in two tabs/users and save both.
Expected result	The system follows an explicit rule: conflict warning, version check or clearly defined last-write behaviour. Silent destructive overwrites must be understood and documented.
Required evidence	Recording plus written explanation of implemented behaviour.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

ADM-ERR-06 - Console and unhandled errors

Applies to	All tested workflows
Team test steps	Run the internal QA suite with browser console open.

Expected result	No uncaught errors, repeated critical warnings or failed requests remain in PASS workflows.
Required evidence	Console screenshot/log attached to the QA package.
Team response	Required: PASS / FAIL / NOT APPLICABLE. If FAIL, explain the cause, owner and target date.

12. Minimum Module Coverage

The team should apply the test cases above to every module included in the pilot. The matrix below defines the minimum coverage. Mark excluded modules as NOT APPLICABLE with an agreed reason.

Module	Navigation/detail	List controls	CRUD/safety	Relationships / special checks
Customers / Users	Required	Search, sort, filter, pagination	Edit/status; delete if supported	Privacy/access; reference record no. 14 flow
Clinics	Required	Required	Create/edit/status/verify	Clinic-Doctor and Clinic-Treatment relationships
Doctors / Experts	Required	Required	Create/edit/status/verify	Clinic relationship; role visibility
Appointments	Required if pilot scope	Required	Edit/cancel/delete/status	Clinic/doctor/customer consistency
Treatments / Services	Required	Required	Full create/edit/publish/deactivate	Devices, tags, translation, app visibility
Devices	Required	Required	Create/edit/deactivate; delete safety	Treatment relation function; protected names
Concerns	Required	Required	Create/edit/delete safety	Search/recommendation impact where applicable
Benefits / Likewise Terms	Required if editable in pilot	Search/list controls as available	Create/edit/delete safety	Impact on tags/search/recommendations
Bookings / Orders / Payments	Required only if pilot scope	Required	Status/error/retry/refund as applicable	Cross-entity integrity and auditability
Support / Help	Required if pilot scope	Required	Create/update/close/escalate	Privacy, ownership and status history
Clinic Panel	Required	Required	Role-specific actions	Only authorised clinic data visible
Doctor Panel	Required	Required	Role-specific actions	Only authorised doctor data visible
Products	NOT APPLICABLE	NOT APPLICABLE	NOT APPLICABLE	Outside current ZYNQ scope

13. Required Evidence Package

- **Build reference:** environment URL, branch, commit/PR and delivery date.
- **Test result matrix:** one row per Test ID with PASS / FAIL / NOT APPLICABLE.
- **Screen recordings:** continuous recordings of the stated workflow; do not provide only static end-state screenshots for navigation tests.
- **Screenshots:** suitable for error, empty, confirmation and language states.
- **Known limitations:** anything mocked, partial, excluded, data-dependent or planned later.
- **Failure plan:** cause, responsible person, estimated effort and target date for every FAIL.

Team response template

Test ID	Status	Environment/ build	Comment	Owner	Target date	Evidence	Remaining limitation
ADM-NAV-01	PASS / FAIL / N/A	Staging build ...	...	...	...	Video/link	...

14. Acceptance Gate

- **Internal QA complete:** All applicable Required tests are executed and evidenced by the development team.
- **No unresolved critical operational failure:** Navigation/state, detail loading, save integrity, permissions and destructive-action safety must pass for pilot-critical modules.
- **Scope is explicit:** Every excluded or partial module is documented and agreed before acceptance.
- **Independent acceptance:** Anna validates representative PASS workflows after the team submits its evidence package.
- **Pilot decision:** Controlled pilot readiness can be confirmed only after the internal QA report and independent acceptance are aligned.

Development delivers and verifies. ZYNQ independently accepts.